

Estonian

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
English


AquaFree Yellow Aspheric Hydrophobic


Hydrophobic Yellow Acrylic Posterior Chamber Intracuclear Lens with UV absorber

Rumex AquaFree Yellow is a foldable posterior chamber hydrophobic acrylic intraocular lens, UV-absorbing optical inner lens used for the replacement of the human crystalline lens in the visual correction of aphakia in the adult patients. Rumex AquaFree Yellow contains blue light filtering chromophore that filters light in a manner the young human crystalline lens does in the 400-475 nm blue light wavelength range.

Rumex AquaFree Yellow is one-piece acrylic lens with square edge, a 6.0 mm biconvex optic, and an overall length of 13.0 mm, manufactured at the following diopter power range: 0 to +10 in 1 diopter increments, and +10 to +30 diopeters in 0.5 diopter increments. Refractive index of the material is 1.5.

Material	Characteristics	UV Cutoff at 10% T
Yellow	UV-Absorber + blue light filter	390 nm

MODE OF ACTION

Rumex AquaFree Yellow posterior chamber intraocular lens is intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This position allows the lens to function as a refractive medium in the correction of aphakia. The effectiveness of the lens in reducing the incidence of retinal disorders has not been established.

INDICATIONS

Rumex AquaFree Yellow posterior chamber intraocular lens is intended for replacement of the human crystalline lens to achieve visual correction of aphakia in adult patients when extracapsular cataract extraction or phacemulsification is performed. The lens is intended for placement in the capsular bag.

PRECAUTIONS

- Do not resterilize the lens by any method.
- Do not store the lens at a temperature greater than 45°C (113°F)
- Use only sterile intraocular irrigating solutions (such as BSS® or BSS PLUS®) to rinse or soak lenses.
- Handle lenses carefully to avoid damage to lens surface or haptics.
- Do not attempt to reshape haptics in any way.
- A high level of surgical skill is required for intraocular lens implantation; the surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on intraocular lens implantation before attempting to implant intraocular lens.

CONTRAINDICATIONS

Patients with any of the following conditions may not be suitable candidates for the intraocular lens implantation because the lens may exacerbate an existing condition, interfere with diagnosis or treatment of a condition, or pose an unreasonable risk to the patient's eyesight. Careful preoperative evaluation and sound clinical judgment should be performed by the surgeon to identify the benefit/risk ratio before implanting a lens in a patient's eye with one or more of the following conditions:

- Choroidal hemorrhage
- Concomitant severe eye disease
- Excessive vitreous loss
- Extremely shallow anterior chamber
- Microphthalmos
- Non-sep related cataract
- Posterior capsular rupture (preventing fixation of IOL)
- Severe corneal dystrophy
- Severe optic atrophy
- Uncontrollable positive pressure
- Zonular separation (preventing fixation of IOL)
- Color vision deficiencies
- Glaucoma
- Chronic uveitis
- Diabetic retinopathy
- Clinically significant macular/RPE changes

WARNINGS

The lens should not be implanted in the following conditions:

- The posterior capsule is ruptured or if a primary capsulotomy is to be performed.
- The peel pouch is found to be damaged or opened.
- Suspected microbial infection.
- Recurrent severe anterior or posterior segment inflammation or uveitis.
- Patients in whom the intraocular lens may affect the ability to observe, diagnose or treat posterior segment disease.
- Surgical difficulties at the time of cataract extraction that might increase the potential for complications (e.g. persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
- A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
- Circumstances that would result in damage to the endothelium during implantation.
- Children under the age of 2 years are not suitable for intraocular lens implantation.

COMPLICATIONS

The following lists the complications which have been associated with the implantation of intraocular lenses (these lists not intended to be all-inclusive)

- Cumulative Adverse Events:
- Hypyema
 - Hypopyon
 - Lens Dislocation
 - Cystoid Macular Edema
 - Pupillary Block
 - Retinal Detachment
 - Intraocular Infection
 - Secondary surgical intervention (excluding retinal detachment and posterior capsulotomy, including, but not limited to the following:
 - a)Iridectomy to reduce a pupillary block
 - b)Vitreous aspiration to reduce a pupillary block
 - c)Repositioning of a lens
 - d)IOL removal to reduce inflammation
 - e)IOL replacement
 - f)Wound leak repair

Persistent Adverse Events:

- Corneal Scar Edema
- Cystoid Macular Edema
- Iritis
- Required IOP requiring treatment

DIRECTIONS FOR USE

- Examine the label on the lens box for proper lens model, diopter power and expiration date.
- Open the lens box to remove the pouched lens and to make sure the lens case information (e.g. power, model and serial numbers) is consistent with the information on the outer box.
- To remove the lens, open the undamaged peel pouch and transfer the case to a sterile environment. Carefully open the case to expose lens. When removing the lens from the case, DO NOT grasp the optic area with forceps. Prior to actual folding process, the lens should be handled by the haptic portion only.
- There are various surgical procedures which can be utilized, and the surgeon should select a procedure which is appropriate for the patient.
- To minimize the occurrence of marks on the lens due to folding, all instruments should be scrupulously cleaned before use.
- DO NOT RESTERILIZE.


HOW SUPPLIED

Rumex AquaFree Yellow is supplied dry, in a lens case packaged in a Tyvek peel pouch and terminally sterilized by ethylene oxide. The lens must be opened only under aseptic conditions (See Instructions For Use above).

EXPIRATION DATE

The packaged Rumex AquaFree Yellow is sterile unless the peel Tyvek seal is damaged or opened. There is a sterility expiration date clearly indicated on the blister pack and the outside box label. Rumex AquaFree Yellow should not be used after that date.

Lithuanian

AquaFree Yellow Aspheric Hydrophobic

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